

=> d his

(FILE 'HOME' ENTERED AT 11:44:25 ON 30 SEP 2004)

FILE 'CAPLUS, USPATFULL, MEDLINE, BIOSIS' ENTERED AT 11:44:57 ON 30 SEP 2004

L1	151433	S	G-PROTEIN?
L2	151433	S	G PROTEIN?
L3	151433	S	L1 (L) L2
L4	61729	S	LUCIFERASE
L5	577	S	BRET
L6	157	S	L5 (L) L4
L7	151433	S	L3 (L) L3
L8	78	S	L3 (L) L6
L9	7745	S	GPCR
L10	51	S	L9 (L) L8
L11	3	S	L10 AND PY <2001
L12	3	S	L8 AND PY<2001

L11 ANSWER 1 OF 3 CAPLUS COPYRIGHT 2004 ACS on STN
TI Detection of β 2-adrenergic receptor dimerization in living cells
using bioluminescence resonance energy transfer (BRET)
PY 2000
AU Angers, Stephane; Salahpour, Ali; Joly, Eric; Hilairret, Sandrine; Chelsky,
Dan; Dennis, Michael; Bouvier, Michel
SO Proceedings of the National Academy of Sciences of the United States of
America (2000), 97(7), 3684-3689
CODEN: PNASA6; ISSN: 0027-8424
AB Heptahelical receptors that interact with heterotrimeric G
proteins represent the largest family of proteins involved in
signal transduction across biol. membranes. Although these receptors
generally were believed to be monomeric entities, a growing body of
evidence suggests that they may form functionally relevant dimers.
However, a definitive demonstration of the existence of G
protein-coupled receptor (GPCR) dimers at the surface of
living cells is still lacking. Here, using bioluminescence resonance
energy transfer (BRET), as a protein-protein interaction assay
in whole cells, we unambiguously demonstrate that the human
 β 2-adrenergic receptor (β 2AR) forms constitutive homodimers when
expressed in HEK-293 cells. Receptor stimulation with the hydrophilic
agonist isoproterenol led to an increase in the transfer of energy between
 β 2AR mols. genetically fused to the BRET donor (Renilla
luciferase) and acceptor (green fluorescent protein), resp.,
indicating that the agonist interacts with receptor dimers at the cell
surface. Inhibition of receptor internalization did not prevent
agonist-promoted BRET, demonstrating that it did not result from
clustering of receptors within endosomes. The notion that receptor dimers
exist at the cell surface was confirmed further by the observation that
BS3, a cell-impermeable crosslinking agent, increased BRET
between β 2AR mols. The selectivity of the constitutive interaction
was documented by demonstrating that no BRET occurred between
the β 2AR and two other unrelated GPCR. In contrast, the
well characterized agonist-dependent interaction between the β 2AR and
the regulatory protein β -arrestin could be monitored by BRET
. Taken together, the data demonstrate that GPCR exist as
functional dimers in vivo and that BRET-based assays can be used
to study both constitutive and hormone-promoted selective protein-protein
interactions.